

TITLE Marlex® HGX-030-01 Mechanical Testing Protocol - Slings

DOCUMENT NUMBER 90740928

PROJECT NAME AND NUMBER Uro/Gyn Sustaining / U0311

AUTHOR/CONTRIBUTORS Carolina Villarreal

APPROVALS Approvals and approval dates are captured through electronic signature within the PDM system.

THESE DOCUMENTS ARE THE PROPERTY OF BOSTON SCIENTIFIC CORP. AND SHALL NOT BE REPRODUCED, DISTRIBUTED, DISCLOSED OR USED FOR MANUFACTURE OR SALE OF APPARATUS WITHOUT THE EXPRESS WRITTEN CONSENT OF BOSTON SCIENTIFIC

Boston Scientific
Marlex® HGX-030-01 Mechanical Testing - Slings
90740928 Rev/Ver. 01
Page 1 of 9

EXHIBIT
34

CONFIDENTIAL
BSCM11500006055

Table of Contents

| | | |
|-----|--|---|
| 1.0 | BACKGROUND | 3 |
| 2.0 | PURPOSE..... | 3 |
| 3.0 | SCOPE..... | 3 |
| 4.0 | REFERENCE DOCUMENTS | 4 |
| 5.0 | BUILD CONDITIONS..... | 5 |
| 6.0 | DATA ANALYSIS AND RECORD RETENTION | 9 |

1.0 BACKGROUND

Phillips Sumika has discontinued production of Marlex® HGX-030-01, a Polypropylene Homopolymer grade. Therefore, Boston Scientific (BSC) obtained the same polymer from a second source (material distributor). While the material has been supplied in the correct Phillips Sumika bags with an identifiable lot number, there was no CoA supplied with this lot, and the distributor is not able to reproduce the CoA. The second lot of material was obtained from the material distributor Emai Plastic Raw Material (Dongguan) Inc

This document outlines the strategy for the mechanical testing that BSC will conduct to ensure the new lot of Marlex® HGX-030-01, received from Emai Plastic Raw Material, is equivalent to the existing lot of Marlex received from Channel Prime Alliance through the current manufacturing process.

2.0 PURPOSE

The purpose of this document is to outline the testing strategy and acceptance criteria that will be used for the lot 6120105 base resin, Marlex® HGX-030-01, to determine lot equivalence when evaluated with the current manufacturing processes at Luxilon, Proxy, AIM and MedVenture Technology Corporation (MTC) to be used in the urethral slings family products.

3.0 SCOPE

The scope of this protocol is to outline the functional material testing that will be conducted for the base resin Marlex® HGX-030-01 received from Emai Plastic Raw Material (Dongguan) Inc at Luxilon, Proxy, AIM, MTC and BSC. The material information is listed in Table 1:

Table 1. Base resin - Material information

| MATERIAL NUMBER | DESCRIPTION | LOT NUMBER | QUANTITY (KG) |
|-----------------|--|------------|---------------|
| 90207538-01 | Polypropylene Resin (Marlex® HGX-030-01) | 6120105 | Natural: 25kg |

All product produced during this build, will be considered non-saleable until final acceptance is documented.

Testing will apply in the representative products for the nominal conditions detailed in Table 2.

Table 2. List of the Representative Product Family

| PROXY / LUXILON / AIM COMPONENT INFORMATION | | |
|--|---------------------------|---|
| BSC NUMBER | PART NUMBER | DESCRIPTION |
| 34106700 | PPS50046 | Polypropylene Monofilament Yarn - 6 mil (Nat) |
| 34114200-02 | 400010008 ES012 - 002 | Advantage Mesh Sheet |
| 90409292-01 | 90409292 | Carrier Base |
| 90409293-01 | 90409293 | Carrier cover |
| MEDVENTURE PRODUCT INFORMATION | | |
| CATALOG NUMBER | UPN PRODUCT NUMBER | DESCRIPTION |
| 850200 | M0068502000 | Advantage™ System |
| 85020005 | M006850200051 | Advantage™ System, Box 5 |
| 850211 | M0068502110 | Advantage Fit™ System |
| 850211 | M0068502111 | Advantage Fit™ System, Box 5 |
| 68504000 | M0068504000 | Obtryx™ System – Curved - Single |
| 68504001 | M0068504001 | Obtryx™ System- Curved - 5-Pack |
| 68505000 | M0068505000 | Obtryx™ System – Halo -Single |
| 68505001 | M0068505001 | Obtryx™ System – Halo - 5-Pack |
| 850-700 | M0068507000 | Solyx™ SIS System – Single |
| 850-700 | M0068507001 | Solyx™ SIS System – 5 Pack |
| 68503000 | M0068503000 | Lynx Suprapubic Sling System |
| 68503001 | M0068503001 | Lynx System 5-Pack |
| 68506000 | M0068506000 | Pre-pubic Sling - Single |
| 68506001 | M0068506001 | Pre-pubic Sling - 5-Pack |

Since all base resin Marlex® HGX-030-01 extrusion parameters will remain the same, no worst case has been selected at Luxilon.

4.0 REFERENCE DOCUMENTS

| Document Number | Description |
|------------------------|--|
| 90036920 | Corporate SOP Capability of Variable Data |
| 90030419 | Corp SOP Design Verification Sampling/ Process Validation Sampling Plans |
| S842767-00 | Corp SOP Risk Analysis |
| S842769-00 | CORP SOP DESIGN VERIFICATION |
| 90356869 | Solyx Risk Analysis Report |
| 90301643 | Solyx Product Specification |
| 90082464 | Lynx Risk Analysis Workbook |
| 90078003 | Lynx System Product Specification |
| 90055034 | Advantage Product Specification |

| Document Number | Description |
|-----------------|----------------------------------|
| 90055035 | Advantage Risk Analysis Workbook |
| 90106687 | Obtryx Risk Analysis Workbook |
| 90104642 | TTFS Product Specification |
| 90141370 | PPS Product Specification |
| 90141908 | Pre-Pubic Risk Analysis Workbook |
| 90732450 | BCA - Marlex HGX-030-01 Mesh |

5.0 BUILD CONDITIONS

Boston Scientific (BSC)

BSC will supply the base resin Marlex® HGX-030-01 to Luxilon and AIM to manufacture the components needed, to perform this protocol.

Luxilon Industries NV (Luxilon)

Luxilon will extrude the base resin Marlex® HGX-030-01 utilizing current approved manufacturing process standards as per their Quality System. Luxilon will then ship the fiber to Proxy for the Mesh processing.

For the purpose of this protocol, a copy of the results obtained will be sent to BSC and the Certificate of Compliance will be sent to Proxy.

PPS50046 - Polypropylene Monofilament Yarn - 6 mil (Nat)

| Characteristic | Unit | Value | Tolerance | Measurement Standard | Key Performance / Process Characteristic |
|--|------|-------|-----------|--|--|
| Denier (Weight in grams per 9,000 meters) | Den | 152 | ± 10 | Spool Average | Key Performance |
| Dtex (Weight in grams per 10,000 meters) | Dtex | 170 | ± 11 | Spool Average | Key Performance |
| Outer Diameter | µm | 152 | ± 16 | WIQA Spool Average | Key Performance |
| Tensile Strength | N | 9.6 | Min | Spool Average | Key Performance |
| Elongation At Break | % | 22 | ± 7 | WIQA Spool Average | Process Characteristic |
| Knot Strength ¹ | N | 7 | ± 1 | WIQA Batch Average | Process Characteristic |
| Shrinkage ¹ | % | 8 | ± 3 | Shrinkage will be measured using a hot air oven at 132°C (270°F) for 10 minutes. Batch Average | Process Characteristic |
| GDP ¹ (Gram Denier Point) | % | 3 | ± 1 | WIQA Batch Average | Process Characteristic |

¹ These are specific process characteristics that are evaluated at Luxilon per their Quality System

Proxy BioMedical (Proxy)

Proxy will receive the fiber from Luxilon, weave and knit the mesh and build the mesh sheets. The mesh sheet (34114200-02/AG) will be produced in accordance with the current manufacturing process standards as per their current Quality System. Traceability will be maintained through lot number assignment and recorded for all components.

All testing will be performed as per test specimen requirement detailed in the applicable test method listed in the table below and mesh sheets will be shipped directly to MTC.

ES012 – 02 Engineering Design Specification Proxy Biomedical Knit-BSC Advantage Mesh

| Test | Min Performance Required | Sample Size & Acceptance Criteria | Specification | Test Method |
|---|---------------------------------|--|--|--------------------|
| Thickness ¹ | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | 0.0026" ± 0.006" | MPS012 |
| Courses per inch (CPI) ² | 95%/ 95% | n = 30, Ppk ≥ 1.21, Pp ≥ 1.24 | 25 -26 | QAP068 |
| Wales per inch (WPI) ² | 95%/ 95% | n = 30, Ppk ≥ 1.21, Pp ≥ 1.24 | 12 ± 1 | QAP068 |
| Aerial Density ³ | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | 100 g/m ² ± 38% | QAP047 |
| Tensile Strength Machine Direction ³ | 95%/ 95% | n = 30, Ppk ≥ 0.74 | ≥37 (lbf) | QAP059 |
| Elongation Test Machine Direction ³ | 95%/ 95% | n = 30, Ppk ≥ 0.74 | 4.1%-14.7%. | QAP060 |
| Device Stiffness ³ | 95%/ 95% | n = 30, Ppk ≥ 0.74 | ≤ 10 (bending length in cm; average face and back) @ 1.1cm width | QAP061 |
| Burst Strength ³ | 95%/ 95% | n = 30, Ppk ≥ 0.74 | ≥ 66 (psi) | QAP062 |

¹ n=30, Ppk ≥ 0.75, Pp ≥ 0.82 per BSC's internal requirements of 95% confidence of 95% reliability (LTPD 5%) for a risk level I characteristic (2-sided spec.) per Corp SOP Design Verification/Process Validation Sampling Plans 90030419.

² n=30 Ppk ≥ 1.21, Pp ≥ 1.24 per BSC's internal requirements of 95% confidence of 95% reliability (LTPD 0.3%) for a risk level II characteristic (2-sided spec.) per Corp SOP Design Verification/Process Validation Sampling Plans 90030419.

³ n=30, Ppk ≥ 0.74 is BSC's internal requirements of 95% confidence of 95% reliability (LTPD 5%) for a risk level I characteristic (1-sided spec.) per Boston Scientific Corp SOP Design Verification/Process Validation Sampling Plans 90030419 & Boston Scientific Corp SOP Design Verification S842769-00. Some of these characteristics may not be RI=I in BSC's dFMEAs listed in section 4.0, but will be tested to those requirements for this protocol.

For the purpose of this protocol, a copy of the results obtained will be sent to BSC and the Certificate of Compliance will be sent to MTC.

AIM Plastics (AIM)

Base resin Marlex® HGX-030-01 will be injected in the mold at AIM during 2 hours, utilizing current approved manufacturing process standards as per their Quality System. Carrier base (90409292-01) and carrier cover (90409293-01) will then be shipped to MTC for processing. Samples will be randomly pulled from the production run from each cavity to perform the following measurements:

CARRIER BASE (90409292-01)

| Characteristic | Min Performance Required | Sample Size & Acceptance Criteria | Specification | Measurement standard |
|-----------------------|---------------------------------|--|----------------------|-----------------------------|
| Critical dimension Q1 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | .040 ± .003 | Vision system |
| Critical dimension Q2 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | .043 ± .003 | Vision system |
| Critical dimension Q3 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | .047 ± .002 | Vision system |
| Critical dimension Q4 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | 2 X .051 ± .005 | Vision system |
| Critical dimension Q5 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | 2 X .058 ± .005 | Vision system |
| Critical dimension Q6 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | .125 ± .005 | Vision system |
| Critical dimension Q7 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | 4 X .030 ± .003 | Vision system |

n=30, Ppk ≥ 0.75, Pp ≥ 0.82 per BSC's internal requirements of 95% confidence of 95% reliability (LTPD 5%) for a risk level I characteristic (2-sided spec.) per Corp SOP Design Verification/Process Validation Sampling Plans 90030419.

CARRIER COVER (90409293-01)

| Characteristic | Min Performance Required | Sample Size & Acceptance Criteria | Specification | Measurement standard |
|-----------------------|---------------------------------|--|----------------------|-----------------------------|
| Critical dimension Q1 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | 4 X .034 ± .002 | Vision system |
| Critical dimension Q2 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | .041 ± .003 | Vision system |
| Critical dimension Q3 | 95%/ 95% | n = 30, Ppk ≥ 0.75, Pp ≥ 0.82 | .104 ± .005 | Vision system |

n=30, Ppk ≥ 0.75, Pp ≥ 0.82 per BSC's internal requirements of 95% confidence of 95% reliability (LTPD 5%) for a risk level I characteristic (2-sided spec.) per Corp SOP Design Verification/Process Validation Sampling Plans 90030419.

For the purpose of this protocol, a copy of the results obtained will be sent to BSC and the Certificate of Compliance will be sent to MTC.

A summary for the general process for the proposed Marlex® HGX-030-01 testing in this protocol is shown in flow chart on Figure 2




Figure 2. Flowchart for Marlex® HGX-030-01 to be used in the urethral sling family products

6.0 DATA ANALYSIS AND RECORD RETENTION

BSC will review all documentation and data generated from each vendor, compile the information, perform data analysis and update this technical report to include all results obtained.



Defining tomorrow, today.™

APPROVAL FOR OFFICIAL DOCUMENT

Document Number:90740928

Version: 01

Lifecycle State: PRE-RELEASED

Effective Date: 18Jan2012

Date Printed: 18Jan2012 16:00 EST

Change Request Number: 1104049

This document is Reference Only unless its version is verified to be current by one of the following methods:

- **Direct lookup in PDM**
- **Shop Floor Paperwork**
- **Officially stamped or otherwise identified as current**

Approval Signature

| Name | Role | Vote | Activity | Date Signed |
|---------------------|-------------|-------------|-----------------|-----------------------|
| Carolina Villarreal | AUTHOR | Accept | Final Review | 18 Jan 2012 15:29 EST |